



CFD-Driven Design of a Low Airflow, Rapid Recovery Containment Control System™




Alexy Kolesnikov
Ray Ryan
Douglas B. Walters

Labs for the 21st century EPA conference. October, 2002
Durham, NC

Efficient Operation. Design Goals.

Task specific custom design loop

Customer requirements



Prototype testing

Airflow modeling and optimization

Efficiently operating enclosure must:

- Maintain high level of containment protection
- Provide a steady balance reading
- Ensure that materials inside the enclosure are undisturbed by airflow
- Provide ergonomic design and ensure ease of access
- Address energy efficiency concerns



Computational Fluid Dynamics (CFD) Principles

Conservation Law System:

Conservation of mass:

$$\frac{1}{\rho} \frac{\partial \rho}{\partial t} + \frac{u_j}{\rho} \frac{\partial \rho}{\partial x_j} + \frac{\partial u_j}{\partial x_j} = 0$$

Conservation of momentum:

$$\rho \frac{\partial u_i}{\partial t} + \rho u_j \frac{\partial u_i}{\partial x_j} - \frac{\partial \sigma_{ij}}{\partial x_j} - \rho b_i = 0$$

Conservation of energy:

$$\rho \frac{\partial e}{\partial t} + \rho u_j \frac{\partial e}{\partial x_j} - \sigma_{ij} \varepsilon_{ij} - \rho s + \frac{\partial q_j}{\partial x_j} = 0$$

σ – Cauchy stress tensor, e – internal energy, s – distributed heat generation,
 ε - strain-rate tensor, q – heat transfer rate by diffusion.


$$\varepsilon_{ij} = \frac{1}{2} \left(\frac{\partial u_i}{\partial x_j} + \frac{\partial u_j}{\partial x_i} \right)$$

- **k- ϵ (wall functions) RANS**
- **eddy viscosity RANS**
- **incompressible flow**
- **time iterated to steady-state**
- **3C**
- **v BC**
- **element study**


Airflow Distribution Parameters

- Enclosure geometry
- Sash height
- Face velocity
- Operator presence
- Room air currents
- Equipment arrangement inside enclosure

Typical fume hood





Typical Flow Sciences enclosure




CFD Airflow modeling



Typical fume hood



Sash door airflow detail



Typical Flow Sciences enclosure




Variable Air Volume Control System

Traditional system design

- Slow or NO response to upset at hood opening
- Indirect flow control measurement (velocity or sash position)
- Primary design focus on energy
- Hood containment not improved by control system performance


Rapid recovery containment control system

- Virtually instantaneous response
- Direct flow control measurement (pressure)
- Primary design focus on containment / safety
- Additional benefit of reduced energy consumption
- Direct coupling between hood design optimization and containment control system performance




CFD Airflow modeling details



Exhaust slots configuration



Airfoil design optimization



Equipment arrangement




Performance Characteristics

Typical fume hood




*Rapid reverse of
airflow at face of the
hood causes loss of
containment while
increasing operator
exposure*



**Flow Sciences, Inc.
Containment System (VBSE_{TM})**



*Engineered enclosure design
and "airfoil" positioning
test to containment beyond
1mg ensuring maximum
operator safety and balance
stability*



Performance Validation

ASHRAE – 110 Tests

<i>Hood ID</i>	<i>Tracer Gas Left Avg. ppm²</i>	<i>Tracer Gas Center Avg. ppm²</i>	<i>Tracer Gas Right Avg. ppm²</i>
<i>2015-1</i>	<0.01	N/A	<0.01
<i>2020-1</i>	<0.01	<0.01	<0.01
<i>2350AC-IP-1</i>	<0.01	N/A	<0.01
<i>2400-1-1</i>	<0.01	<0.01	0.01
<i>3100-1</i>	<0.01	N/A	<0.01
<i>8500-1</i>	<0.01	<0.01	<0.01

Standards



ASHRAE 0.1 ppm ~600,000ng/m³

PEL < 1000 ng/m³



Flow Sciences, Inc. < 5ng/m³

Potential Design Solutions

Back baffle design variation




Inflow redirection solution





Final Design Detail



Velocity magnitude profile



Positive horizontal velocity component






*Final design
Patent pending*



ASHRAE 110 test results at 40 ft/m

Escape characteristics




<i>Hood ID</i>	<i>Tracer Gas Left Average ppm</i>	<i>Tracer Gas Right Average ppm</i>
<i>New design</i>	0.00	0.00


**ASHRAE
standard 0.1 ppm!!**

ASHRAE 110 test results at 80 ft/m


Escape characteristics



Left



Right







<i>Hood ID</i>	<i>Tracer Gas Left Average ppm</i>	<i>Tracer Gas Right Average ppm</i>
<i>New design</i>	0.00	0.00



**ASHRAE
standard 0.1 ppm!!**

300 ft/m metal plate Blow-by Tracer Gas Test

40 ft/m




80 ft/m




Additional Tracer Gas Challenge

Work area concentration validation



Test
position




	<i>New Design ppm average</i>	<i>Traditional Design ppm average</i>
<i>40 ft/m</i>	3	>120
<i>80 ft/m</i>	3	60

Significant
> 95 – 97.5% reduction of “roll”

Occupied enclosure simulation

- Dynamic equipment re-arrangement test
- Dynamic hood operation test



Engineered for Safety and Performance

Rapid Recovery VAV System Validation

Computer Driven
Air Flow Optimization

Industry Standard
Compliance

