



CFD-Driven Optimization of Room Airflow Patterns And Their Effect on Hood Containment for Flexible Laboratories of the Future



Alexy Kolesnikov
Ray Ryan
Douglas B. Walters

Labs for the 21st century EPA conference. October, 2003
Denver, CO

Modern Laboratory Containment Demands

Balance enclosure



Powder handling



Robotic/High throughput technologies



Containment Performance Characteristics

Work area airflow distribution

- Ensure non-turbulent, unidirectional airflow distribution
- Provide ergonomic design and ensure ease of access
- Address energy efficiency concerns

Operator presence, laboratory operation

- Instantaneous disturbances at the face of the hood
- Dynamic airflow distribution in the Laboratory


Laboratory airflow patterns

- Hood position
- Diffuser blanking
- Diffuser/hood position
- Diffuser/hood separation
- Transfer grilles
- Make up air distribution




CFD Airflow modeling




Typical fume hood



Typical Flow Sciences enclosure




Sash door airflow detail




CFD Airflow modeling details



Exhaust slots configuration



Airfoil design optimization




Equipment arrangement




Variable Air Volume Control System

Rapid recovery containment control system

- Virtually instantaneous response
- Direct flow control measurement (pressure)
- Primary design focus on containment / safety
- Additional benefit of reduced energy consumption
- Direct coupling between hood design optimization and containment control system performance




Prototype development



Final design

- *flexible*
- *portable*
- *efficient*

Patent pending



- Recirculation-free laminar airflow for maximum safety
- Variable Air Volume control system for sash door cross-current compensation
- Low airflow energy efficient design
- Ergonomic “telescoping” sash configuration



Room airflow modeling and hood containment optimization

Containment performance parameters

- Room geometry
- Mechanical HVAC equipment
- Diffuser type and placement
- Operational procedures





Room airflow modeling and hood containment optimization

Current Laboratory guidelines are designed to promote containment **inside** the hood



They **DO NOT** address the fact the containment performance is linked with the air movement in the laboratory around the sash opening



The problem is compounded by the modern Laboratory **flexibility** demands




Experimental tests are costly, provide limited data and identify the problem only **after** the laboratory has been designed and built

Room airflow modeling and hood containment optimization


Possible recommendations

- Hood position
- Diffuser blanking
- Diffuser/hood position
- Diffuser/hood separation
- Transfer grilles
- Make up air distribution
- Hood separation, same wall
- Hood separation, opposite walls
- Hood separation, perpendicular walls




Room airflow. Face velocity distribution

Velocity vector



$$\bar{\mathbf{V}} = \bar{\mathbf{V}}_x + \bar{\mathbf{V}}_y + \bar{\mathbf{V}}_z$$

$$V^2 = V_x^2 + V_y^2 + V_z^2$$



V_z vertical cross flow

V_y transverse cross flow

V_x (measured face velocity)

Ideal airflow distribution




$$\mathbf{V} = \mathbf{V}_x$$

$$\mathbf{V}_y = 0$$

$$\mathbf{V}_z = 0$$


Model Configurations

- Hood in isolation (ideal testing conditions)
- Two hoods, opposite walls
- Two hoods, adjacent walls




Single Hood in Isolation


Velocity magnitude. Top view




Velocity magnitude. Cross view



Inflow distribution




Cross flow distribution




Hood separation. Opposite walls


Velocity magnitude. Top view




Velocity magnitude. Cross view



Inflow distribution, Hood 1




Cross flow distribution, Hood 1




Hood separation. Opposite walls


Velocity magnitude. Top view




Velocity magnitude. Cross view



Inflow distribution, Hood 2




Cross flow distribution, Hood 2




Hood separation. Perpendicular walls

Velocity magnitude. Top view




Velocity magnitude. Cross view



Inflow distribution, Hood 1


Normalized Magnitude



Node number

Cross flow distribution, Hood 1


Normalized Magnitude




Node number

Hood separation. Perpendicular walls

Velocity magnitude. Top view




Velocity magnitude. Cross view



Inflow distribution, Hood 2


Normalized Magnitude



Node number

Cross flow distribution, Hood 2

Normalized Magnitude



Node number



Potential recommendations. Additional considerations

Hood positioning recommendations

- Avoid positioning hoods next to a perpendicular wall
- Avoid placing hoods close to each other at perpendicular walls
- Provide sufficient space at the hood face opening to ensure positive incoming flow distribution
- Avoid positioning hoods next to doors, windows and high traffic areas

Additional considerations

- Diffuser blanking
- Diffuser/hood position
- Diffuser/hood separation
- Transfer grilles
- Make up air distribution



CFD room airflow modeling

CFD simulation protocols are very detailed

- Mathematics
- Physics
- Computer science

CFD provides data otherwise absolutely unavailable

- Repeatability
- Parametric variation
(hood, furniture, supply, exhaust)
- Density of data





Room airflow modeling and hood containment optimization

- Safety**


Predict and optimize the containment performance of the complete laboratory/hood airflow system

- Flexibility**

Be able to make changes to the entire system layout knowing the exact effect of such variations

- Energy conservation**

Use only enough air changes per hour to satisfy the requirements. Hood face velocity optimization



↓

Project Cost Benefits

↓

- Design optimization prior to construction or mock-up
- Project time reduction via established CFD data management protocols utilization
- Energy demand minimization

